
## **ELEVATE JUMP**™

## ELEVATE





Incline Release Lever

to move between

seven levels of resistance.

**JumpOMeter** allows quick visual measurement progression of plyometrics.

Variable Band

allows for an

Resistance (VBR)

additional 10-70 lbs of resistance.

rement ometrics.

Curved Platform follows anatomic curve of hip to foot ratio ensuring proper biomechanics throughout plyometric movements.

Slide Distance Regulator

## to control range of motion.



Anchor Feet with Bolt Holes for safety and stability.

**Safety Lock** enables users begin exercise in a safe and comfortable position.

## PRODUCT SPECIFICATIONS

IN USE 104" x 33" x 68" [L/W/H] (2.6 m X .8 m X 1.7 m)

**FLOOR SPACE** 24 square feet

UNIT WEIGHT 230 lbs [105 kg]

**CONSTRUCTION** Steel with steel reinforced

extruded aluminum rails

(2.2 square meters)

**INSTALLATION** Designed to be bolted to the floor

the moor

**CAPACITY** Maximum user weight capacity

of 400 lbs [180 kg]

BODYWEIGHT 50% RESISTANCE

50% - 80% BWR

**BANDS** Additional 10 – 70 lbs VBR

WARRANTY Frame: 5 years
Upholstery: 90 days

Rubber & Moving Parts: 1 year

The ELEVATE Jump™ builds explosive power and develops lower body muscle mass with two primary movements, the squat and the squat jump. Uniquely designed to facilitate plyometric exercises leveraging both linear bodyweight resistance and band variable resistance, the Jump produces both concentric and eccentric loading of the muscles − allowing users to accelerate and decelerate quickly during a squat jump in a safe, controlled manner. Able to accommodate all levels of fitness, the ELEVATE Jump supports the spine during traditional squat exercises so that users can comfortably control the descending and ascending phases of the squat, while providing an opportunity for advanced variations including single-leg squats and staggered stance squats.

- **1.** Provides concentric and eccentric force utilizing both linear bodyweight resistance (BWR) and variable band resistance (VBR).
- **2.** Four resistance bands can be engaged to allow for an additional 10-70 lbs resistance in 10 lb increments.
- **3.** Pneumatic brake protects joints during plyometric movements.
- **4.** Angled glideboard provides comfort and support for the upper body.
- **5.** Includes link to exercise library with additional and advanced exercise options.